

### Lesson 3a – dilution before titration-solutions

The acetic acid (ethanoic acid) concentration of a brand of vinegar is to be determined using volumetric analysis. A 20.00 mL aliquot is taken from the original bottle of vinegar and placed in a 200 mL volumetric flask and made to the mark using distilled water.

A volume of 25.00 mL was transferred from the volumetric flask to a 100mL conical flask and titrated to the end point using a standard solution of 0.201 M  $\text{Na}_2\text{CO}_3$ . An average titre of 20.16 mL was obtained. Find the concentration of acetic acid in the original sample in %m/v.



- a) Write the balanced overall equation for the reaction taking place in the conical flask between the ethanoic acid and the  $\text{NaHCO}_3$ .



- b) Find the mol of  $\text{Na}_2\text{CO}_3$  in the average titre

$$\text{Mol of } \text{Na}_2\text{CO}_3 = C \times V = 0.201 \text{ mol/L} \times 0.02016 \text{ L} = 4.05 \times 10^{-3} \text{ mol}$$

- c) Find the mol of acetic acid in the conical flask.

$$\text{Mol of acetic acid} = 2 \times \text{mol of } \text{Na}_2\text{CO}_3 = 8.10 \times 10^{-3} \text{ mol}$$

- d) Find the concentration, in mol/L, of acetic acid in the volumetric flask

*Since the  $8.10 \times 10^{-3} \text{ mol}$  of acetic acid came from a 25.00 mL sample of the solution in the volumetric flask, we can calculate the concentration of acetic acid in the 25.00 mL sample.*

$$\Rightarrow 8.10 \times 10^{-3} \text{ mol} / 0.0250 = 0.324 \text{ M}$$

*concentration of acetic acid in the volumetric flask = concentration of acetic acid in the 25 mL sample*

$$\text{hence concentration of acetic acid in the volumetric flask} = 0.324 \text{ M}$$

- e) Find the concentration in mol/L in the original undiluted sample

*Since the original sample was diluted we can use the formula*

$$C_1 V_1 = C_2 V_2$$

$$C_1 = C_2 V_2 / V_1$$

$$\Rightarrow C_1 = 0.324 \times 0.200 / 0.0200 = 3.24 \text{ M}$$

- f) Find the concentration of acetic acid, in %m/v, in the original sample .

*Convert mol/L to %m/v*

$$\Rightarrow 3.24 \text{ mol/L} = (3.24 \times 60.0 \text{ g}) / 1000 \times 100 = 19.44\% \text{ m/v}$$

Unpack the information by drawing a flow diagram.

Convert mol/L to %m/v  
19.44% m/v

20.00 mL  
original  
sample

200 mL volumetric  
flask



25.00 mL  
aliquot is  
taken



An average titre  
20.16 mL of a  
0.201 M  $\text{Na}_2\text{CO}_3$ .



Concentration of  $\text{CH}_3\text{COOH}$  in the  
volumetric flask is 0.324 M  
This was diluted from the original  
sample in a ratio 1:10.

$$C_1V_1 = C_2V_2$$

$$C_1 = C_2V_2/V_1 = 0.324 \cdot 0.200 / 0.0200$$

$$C_1 = 3.24 \text{ M}$$

$8.10 \times 10^{-3}$  mol of  $\text{CH}_3\text{COOH}$  in the  
conical flask that came from 25.00 mL  
of solution from the volumetric flask  
Concentration of  $\text{CH}_3\text{COOH}$   
 $= 8.10 \times 10^{-3} \text{ mol} / 0.025 = 0.324 \text{ M}$

$4.05 \times 10^{-3}$  mol of  
 $\text{Na}_2\text{CO}_3$