

Using spectroscopy to determine purity

Worksheet – Use of spectroscopy in qualitative and quantitative analysis.


Steps for verifying purity using IR spectroscopy:

1. **Obtain the reference spectrum:**
 - Obtain the IR spectrum of the pure substance from the literature.
2. **Identify characteristic troughs:**
 - Note the specific absorption troughs and their corresponding wavenumbers in the reference spectrum that are characteristic of the pure substance.
3. **Analyze the sample's IR spectrum:**
 - Obtain the IR spectrum of the sample in question.
4. **Compare Spectra:**
 - Compare the sample's IR spectrum with the reference spectrum.
 - Check if all the characteristic troughs of the pure substance are present in the sample's spectrum at the same wavenumbers.
5. **Look for additional troughs:**
 - Identify any additional troughs in the sample's IR spectrum that are not present in the reference spectrum, use the 2024-2027 data booklet.
 - Unexpected troughs may indicate the presence of impurities.
6. **Interpret the troughs:**
 - Analyze the unexpected troughs to identify possible impurities by comparing them with known IR spectra of other compounds.
7. **Draw Conclusions:**
 - Determine if the sample's spectrum matches the reference spectrum, indicating purity.
 - If there are additional troughs, assess their intensity and nature to estimate the extent and type of impurities.

Example 1.

Analyse the purity of a sample of ethanol.

1. Obtain the reference spectrum:




2. Identify characteristic troughs:

Note the specific absorption troughs and their corresponding wavenumbers in the reference spectrum that are characteristic of the pure substance. Use the 2024-2027 data booklet.

Characteristic ranges for infrared absorption			
Bond	Wave number (cm ⁻¹)	Bond	Wave number (cm ⁻¹)
C=O (amides)	1630–1680	C–H (alkanes, alkenes, arenes)	2850–3090
C=O (aldehydes)	1660–1745	O–H (acids)	2500–3500
C=O (acids)	1680–1740	O–H (alcohols)	3200–3600
C=O (ketones)	1680–1850	N–H (amines and amides)	3300–3500
C=O (esters)	1720–1840		

IR from the 2024-2027 data booklet




O-H bond at 3200-3600

3. Analyze the sample's IR spectrum:

Obtain the IR spectrum of the sample in question.

4. Compare spectra:

All the characteristic signals of ethanol are present in the sample IR.



5. Look for additional troughs and interpret the signals:

An unexpected trough at 1660 indicates the presence of an aldehyde C=O according to the data booklet.

6. Draw conclusions:

Using the information provided and the references from the data booklet it is likely that the sample of ethanol is contaminated with an aldehyde.

Steps for verifying purity using mass spectroscopy:

1. Obtain the reference mass spectrum:

- Find and study the mass spectrum of the pure substance from the literature.

2. Identify key peaks:

- Note the key peaks (m/z values) and their relative abundances in the reference spectrum that are characteristic of the pure substance.

3. Analyse the sample's mass spectrum:

- Obtain the mass spectrum of the sample in question.

4. Compare Spectra:

- Compare the sample's mass spectrum with the reference spectrum.
- Check if all the key peaks of the pure substance are present in the sample's spectrum at the same m/z values.

5. Look for Additional peaks:


- Identify any additional peaks in the sample's mass spectrum that are not present in the reference spectrum.
- Unexpected peaks may indicate the presence of impurities.

6. Interpret the peaks:

- Analyse the unexpected peaks to identify possible impurities by comparing them with known mass spectra of other compounds.
- Look for fragment ions and their patterns to gain more information about possible impurities.


Example 2 – Ethanol was produced via a chemical pathway where the feedstock was ethene. The final product was analysed for impurities and extent of reaction using mass spectroscopy. Analyse the purity of a sample of the ethanol, whose MS is given to you, and comment on whether this pathway yields 100% ethanol. The MS of pure ethanol is also given provided.

1. Obtain the reference mass spectrum:




2. Identify key peaks:

- m/z 46 molecular ion $[\text{CH}_3\text{CH}_2\text{OH}]^+$
- base peak at m/z 31 $[\text{CH}_3\text{OH}]^+$



3. Analyse the sample's mass spectrum:

- peaks at m/z 46 and 31 are still there indicating ethanol is present.
- new peaks at m/z 44
- peak at m/z 15 possibly a $[\text{CH}_3]^+$, has increased in size indicating the presence of a CH_3 in the impurity.
- peak at m/z 28 is now the base peak indicating the presence of a stable fragment in the impurity, such as $[\text{CH}_3\text{CH}]^+$
- m/z 44 now shows a strong signal most likely the parent ion of the impurity such as $[\text{CH}_3\text{CHO}]^+$




4. Possible impurity might be ethanal (CH_3CHO)



Obviously ethanal is the precursor to the formation of ethanol and makes sense that we have identified an unused amount of ethanal still in the final product mix. The reaction does not produce a 100% yield.



1 ↓ X



2 ↓ $\text{Cr}_2\text{O}_7^{2-}/\text{H}^+(\text{aq})$



3 ↓ "Y"




Figure 1 – Reaction pathway of propene.

1. Consider the reaction pathway shown in figure 1. The IR spectra of the different products formed are shown for wave numbers between 1500 cm^{-1} to 4000 cm^{-1} . Samples of each substance A, B, C and D were analysed after each step with IR spectroscopy. Substances A, B, C and D can be mixtures or pure substances.

a. What are the reagents and conditions represented by "X"?

b. Consider substance "A".

i. To what class of compounds does "A" belong to?

ii. Justify your answer with reference to the IR spectrum to question i. above..

iii. Is "A" a pure substance? _____

iv. Justify your answer to question iii. above, with reference to the reaction pathway shown.

c. What are the reagents and conditions represented by "Y"?

d. Name (IUPAC) the compound formed during reaction "3".

e. Consider substance "D".

i. Is substance "D" a pure substance or a mixture?

ii. Justify your answer to question i. above with reference to the IR spectrum of substance "D".

-
- Figure 2 -MS of pure propan-1-ol




Figure 3 – MS of sample

This image shows a single sheet of white paper with horizontal blue or grey ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

3. You are provided with the IR spectrum of an unknown ethanol solution in fig 4. The absorbance (measured as trough height) of the O-H stretch band around 3400 cm^{-1} for the ethanol sample is 0.75 cm . A calibration curve is to be created using solutions of ethanol, and the following transmittance (trough heights) data was obtained, table 1, for the O-H stretch




Figure 4 – IR transmittance of ethanol

Concentration (mg/L)	Trough height at (3400 cm^{-1}) in cm
10	0.2
20	0.4
30	0.6
40	0.8
50	1.0

Table 1 – concentration vs trough height

a. Using the data in table 1 plot a calibration curve on the graph paper provided on the next page and draw the line of best-fit.


b. Using the calibration curve, give the concentration of ethanol in the sample in:

i. mg/L

ii. mol/L

iii. ppm

iv. % v/v (density of ethanol 0.789 g/mL)



4. You are provided with a mass spectrum of an unknown ethanol sample. A calibration curve was previously created using known ethanol concentrations and the absolute, rather than the relative heights, were measured. The following data was obtained using peak height for the m/z 31 ion:





Figure 5 – MS of ethanol absolute values.

Concentration (mg/L)	Absolute value of m/z 45 (mm)
10	0.24
20	0.51
30	0.74
40	1.02
50	1.26

Table 2 – concentration vs absolute value of peak.

- Using the data in table 1 plot a calibration curve on the graph paper provided on the next page and draw the line of best-fit.
- What fragment formed the peak at 45 m/z ? _____
- Using the calibration curve, give the concentration of ethanol in the sample in:
 - mg/L
 - mol/L
 - ppm
 - % v/v (density of ethanol 0.789 g.mL)



5. Acetaldehyde (C_2H_4O) is formed via an oxidation reaction with ethanol. Consider the mass spectra of pure ethanol and a product sample of ethanol.

a. Consider the sample MS. Identify the fragments that formed peaks at the following m/z values from the simple splitting of the parent ion eg just one bond break.

- 15 _____

- 29 _____

- 44 _____

- 45 _____

- 43 _____

- 46 _____

- 31 _____








Figure 6 – mass spectra of pure ethanol and sample.

- b. Has the reaction to convert ethanol into ethanal resulted in a 100% yield? Explain with evidence from the spectra given in figure 6.

- c. Why would the fragment at m/z 15 not be an accurate indicator of contamination?

6. Consider the three MS of three organic molecules (alcohols) with the molecular formula $C_4H_{10}O$. Identify compounds A, B and C. Justify each choice with reference to the MS of each compound by completing table 1 below.



Molecule (IUPAC)	Skeletal structure	Identifying m/z signal/s due to simple, one-bond-break, splitting.	Fragment/s.
A			
B			
C			

Table 1