Lesson 1 IR spectroscopy

Click to revise IR

It is expected that students have revised IR by visiting the link above.

Molecules are not rigid, inflexible structures but can bend and stretch in many different ways, as the image of a water molecule shows on the link above. It is best to view the online IR section so that you can fully appreciate the movement of bonds which causes absorption of IR radiation. Covalent bonds in organic molecules are not rigid and behave more like springs. At room temperature, organic molecules are always changing shape ever so slightly, as their bonds undergo constant stretching, bending and twisting. These complex vibrations are not random and each one occurs at

specific frequencies. This means that a molecule can only stretch and bend at certain 'allowed' frequencies. Take the carbonyl group from ketone, as shown on the right. It's stretching occurs at 5.15 X 10²³ times every second at room temperature. If the molecule is exposed to electromagnetic radiation that matches the frequency of the carbonyl stretch, it will absorb this energy and increase the amount of stretch, but not the frequency of



stretching. That is, each stretch will be longer while the rate remains at of 5.15 X 10²³ per second. This takes the bond to a higher energy state. It turns out that the energy that bonds can absorb fall within the IR region of the electromagnetic spectrum.

The frequency at which a particular bond type absorbs energy depends on the strength of the bond and the mass of each of the atoms on either side of the bond.

Molecules such as Cl₂ are IR inactive because they have no dipoles. Carbon dioxide, for example, although symmetrical, is IR active. When its bonds undergo an asymmetric stretch or bend the dipole moment of the molecule changes at a certain rate per second, known as its frequency measured in Hz. If this changing dipole frequency matches a frequency in the IR electromagnetic region then absorption of that particular IR frequency takes place. You will notice that IR absorption spectra graph absorbance or transmittance, on the Y-axis, against wavenumber, on the X-axis. Wavenumber is given as how many waves can fit into one centimetre hence its units are cm⁻.

Frequency and wave number are related but the explanation is beyond the scope of this course.

An IR spectrum reveals the type of bond present and can also be used to identify a compound by its fingerprint region.



Example 1 Consider the IR spectrum on the right.

A researcher has two unlabelled bottles A and B that contain either pentanoic acid or pentanol. Identify the compound in bottle A given its IR spectrum on the right.



A strong absorption at 1700cm⁻ indicates an acidic carbonyl group. The broad absorption band between 2500 – 3500 indicates an acidic OH.



The substance in bottle B gives the IR spectrum shown on the right.

Identify the key feature that identifies the compound as the alcohol.



An absorption band between 3200 cm⁻ and 3600 cm⁻ is indicative of an alcohol OH group.



 Which of the following molecules are IR active? Give a reason. CO₂, Cl₂, CH₄.

 Cl_2 is the only moleuycle out of the three which IR inactive. In order for a molecule to absorb in the IR region its dipole moment must change as the molecule vibrates or rotates. Using the molecule HCl as an example, the charge distribution betwee hydrogen and chlorine is not evenly distributed since chlorine is more electronegative than hydrogen and so has a higher electron density surrounding it. This creates a large dipole moment and makes the bond polar. As the molecule vibrates, fluctuations in its dipole moment occurs. This changing dipole interacts with the electromagnetic radiation and absorbs in the IR region at a frequency identical to that of the chanign dipole moment.

Molecules such as O_2 , H_2 , F_2 , o, $r N_2$ can not produce a changing dipole moment when their bonds undergo rotational or vibrational changes. Ss a result, they do not absorb IR radiation.

 A chemist wanted to identify the contents of a container. The IR spectrum of a pure sample taken from the bottle is shown on the right.

Which of the two suspected compounds is present in the bottle.

Butanol or butanal. Explain.

Butanal has the characteristic carbonyl absorption band at about 1700 cm⁻





The butanol, on the other hand, should have a characteristic band at 3200-3600 cm⁻ indicating an alcohol OH group.



- 3) The frequency at which a particular bond absorbs IR radiation depends on what factors? *The mass of the atoms on either side of the covalent bond and the strength of the bond itself.*
- 4) The IR spectrum of a compound with the molecular formula C₄H₈O₂ is shown on the right. Identify the molecule. Give reasons.
 Butanoic acid

